AJMB

Avicenna Journal of Medical Biochemistry

Avicenna J Med Biochem, 2022; 10(2):84-89. doi:10.34172/ajmb.2022.2358

http://ajmb.umsha.ac.ir

Angiotensinogen Gene M235T and T174M Polymorphisms in Diabetic Nephropathy in a Bangladeshi Population

Md. Mizanur Rahman^{1,2*}, Shahidul Islam Salim³, Imran Khan⁴, Khondoker Moazzem Hossain², Liaquat Ali^{5,6}, Zahid Hassan⁷

Supplementary Table S1: Reagent 1 Used for Glucose Concentration Measurement Assay

Contents	Initial Concentration	of
	Solution	
Buffer:		
Phosphate Buffer	0.1 mol/I, pH 7.0	
Phenol	11 mol/I	
GOD-PAP Reagent:		
4-aminophenazone	0.77 mmol/I	
Glucose oxidase	≥1.5 kU/I	
Peroxidase	≥1.5 kU/I	
Standard:		
Glucose	5.55 mmol/I (100 mg/dI)	
Additional Reagent:	`	
Uranyl Acetate	0.16 %	

¹Department of Medicine, Division of Nephrology, Faculty of Medicine and Dentistry, University of Alberta, Edmonton, T6G2G3, Canada

²Biotechnology and Genetic Engineering Discipline, Life Science School, Khulna University, Khulna-9208, Bangladesh

³Department of Nephrology, Bangabandhu Sheikh Mujib Medical University, Dhaka-1000, Bangladesh

⁴Incepta Vaccines Limited, Dhaka, Bangladesh

⁵Department of Biochemistry & Cell Biology, Bangladesh University of Health Sciences, 125/1 Darus Salam, Mirpur-1, Dhaka-1216, Bangladesh

⁶Department of Research & Development, Pothikrit Institute of Health Studies, Dhaka, Bangladesh

⁷Department of Physiology & Molecular Biology, Bangladesh University of Health Sciences, 125/1 Darus Salam, Mirpur-1, Dhaka-1216, Bangladesh

Supplementary Table S2: Reagent 2 Components (Consists of Picric Acid and Sodium Hydroxide, mixed in Equal Volume) for Creatinine Concentration Measurement Assay

Contents	Concentration of Solution
Creatinine standard	177 umo1/I (2 mg/dI)
Picric acid surfactant	35 mmol/I
Sodium hydroxide	0.32 mol/I